

Tissue Drying Oven TDO 66

- Quickest drying of large specimen quantities



Increase image

Key product features

Slide drying oven with electronic control and digital display with hot air ventilation system using a fan. The new model of **Medite**[®] slide drying oven **TDO 66** with electronic temperature control, timer and digital display provides unique features. Although designed primarily for histology, because of its flexibility this quick dryer is applicable also in other laboratories.

The particular advantage of the **Medite**[®] **TDO 66** is the extremely short warming up and drying time. Using the thermostat controlled fan, heat is regularly distributed inside the chamber and the set temperature is achieved already in about 1 minute. Total drying process is completed in only 5 to 10 minutes depending on the desired temperature.

Equipment: ready to use with two lattice inserts.

Technical data

Programming:

programmable drying time from 1 min. to 99 h 59
or operation with timer

Dimensions overall (W / D / H):

450 x 300 x 470 mm

Drying chamber (W / D / H):

400 x 270 x 240 mm

Temperature:

+20° C to +80° C

Weight:

17 kg

Power supply:

115-230 V / 50/60 Hz / 1500 VA

[Product information US SDF 2008.pdf](#)



1.3 M

download detailed product information